

Pe Mouse Anti-Human Cd56/Ncam1 Monoclonal Antibody For Flow Cytometry - P13591

Item Number: CM0004433



Introduction

PE Mouse anti-Human CD56/NCAM1 monoclonal antibody for flow cytometry applications. Detects the 95 kDa neural cell adhesion molecule 1 (NCAM1), which is involved in neuron adhesion and viral entry. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD56/NCAM1 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
Gene ID	4684
Immunogen	Recombinant protein of human CD56/NCAM.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag