

Pe Mouse Anti-Human Tcr Vδ2 Monoclonal Antibody For Fc - P04234

Item Number: CM0011013



Introduction

Optimized for flow cytometry, this PE-conjugated mouse monoclonal antibody against human TCR Vδ2 delivers reliable T-cell subset analysis and high specificity for human samples in immunological research, ideal for gamma/delta T cell studies.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human TCR Vδ2 mAb
Remarks/Alias	TCRD; TRD@; TCRDV1
Species	Human
GeneID	6964
Immunogen	Recombinant protein of human TCR Vδ2
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04234
Protein Weight	19kDa
Shipping	Ice bag