


Pe Mouse Anti-Rabbit IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01870

Item Number: CM0006472



Introduction

PE-conjugated mouse monoclonal antibody targeting the Fc region of rabbit IgG. Designed for flow cytometry detection of rabbit immunoglobulins in immunological research and diagnostic assays.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Rabbit IgG (Fc) mAb
Species	Rabbit
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 96-323 of rabbit IgG (Fc) (P01870)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P01870
Protein Weight	35kDa
Shipping	Ice bag