

Pe Rabbit Anti-4-Hydroxynonenal (4-Hne) Monoclonal Antibody For Flow Cytometry

Item Number: CM0000056

Introduction

PE-conjugated rabbit monoclonal antibody targeting 4-Hydroxynonenal (4-HNE), a lipid peroxidation marker. Suitable for intracellular flow cytometry (FC). Species-independent reactivity. Supplied by CamelBio for IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-4-Hydroxynonenal mAb
Remarks/Alias	4-HNE
Immunogen	Chemical compounds corresponding to 4-Hydroxynonenal.
Source	Rabbit
Category	Other Modified Antibodies
Application	FC (intra)
Cross Reactivity	Species independent
Shipping	Ice bag