

Pe Rabbit Anti-Chicken Cd25 Monoclonal Antibody For Fc - Q90Wj4

Item Number: CM0004552

Introduction

PE-conjugated rabbit monoclonal antibody targeting chicken CD25 (IL-2 receptor alpha), validated for flow cytometry. Ideal for avian immunology research, Treg studies, and immune tolerance investigations.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Chicken CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43; p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Mouse
Gene ID (Human)	3559
Gene ID	395294
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-190 of chicken CD25 (NP_989927.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Chicken
SWISS	Q90WJ4
Protein Weight	24 kDa
Shipping	Ice bag