

Pe Rabbit Anti-Dog Cd19 Monoclonal Antibody - Cd19

Item Number: CM0008193

Introduction

High-quality PE-conjugated rabbit monoclonal antibody against canine CD19, specifically validated for flow cytometry. Suitable for B-cell immunophenotyping and immune monitoring in dog models. Optimized for reliable performance in IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Dog CD19 mAb
Species	Dog
Gene ID (Human)	930
Gene ID	607898
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 20-299 of dog CD19 (XP_013969991.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Dog
Shipping	Ice bag