

Pe Rabbit Anti-G4S Linker Mab - G4S Linker

Item Number: CM0009696

Introduction

PE-conjugated rabbit monoclonal antibody targeted to G4S linker peptide. Validated for flow cytometry (FC) with species-independent reactivity, suitable for detecting G4S-tagged proteins in research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-G4S Linker mAb
Immunogen	Synthetic peptide (A synthesized peptide derived from G4S Linker)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag