

Pe Rabbit Anti-Goat IgG (Fc) Fragment Mab - Flow Cytometry Antibody

Item Number: CM0012484

Introduction

PE-conjugated rabbit monoclonal antibody specific to goat IgG (Fc) fragment. Designed for flow cytometry, offering high specificity and minimal cross-reactivity with other species. Ideal for immunoassay detection of goat antibodies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit Anti-Goat IgG (Fc) fragment mAb
Species	Goat
Immunogen	Recombinant protein of goat IgG (Fc) fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Goat
Shipping	Ice bag