


Pe Rabbit Anti-Human C5Ar1/Cd88 Monoclonal Antibody For Fc - P21730

Item Number: CM0009376



Introduction

High-quality PE-conjugated rabbit monoclonal antibody targeting human C5AR1/CD88. Specified for flow cytometry (FC). For inflammation, chemotaxis, and GPCR signaling research. CamelBio IVD raw material supply.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human C5AR1/CD88 mAb
Remarks/Alias	C5A; C5AR; C5R1; CD88
Species	Human
GeneID (Human)	728
GeneID	728
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human C5AR1/CD88 (NP_001727.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21730
Protein Weight	39kDa
Shipping	Ice bag