


Pe Rabbit Anti-Human Calreticulin Monoclonal Antibody For Intracellular Flow Cytometry - P27797

Item Number: CM0008827



Introduction

PE-conjugated rabbit anti-human Calreticulin monoclonal antibody for intracellular flow cytometry (FC). Targets the ER chaperone CALR (P27797), involved in protein folding, calcium binding, and quality control. Reacts with human.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Calreticulin mAb
Remarks/Alias	RO; CRT; SSA; cC1qR; HEL-S-99n
Species	Human
GeneID (Human)	811
GeneID	811
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 335-417 of human Calreticulin (NP_004334.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P27797
Protein Weight	48kDa
Shipping	Ice bag