

Pe Rabbit Anti-Human Casr Monoclonal Antibody For Fc - P41180

Item Number: CM0008695



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CaSR (extracellular calcium-sensing receptor; SWISS: P41180, ~121 kDa). Validated for flow cytometry; supplied as IVD raw material for diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CaSR mAb
Remarks/Alias	CAR; FHH; FIH; HHC; EIG8; HHC1; NSHPT; PCAR1; hCasR; GPRC2A; HYPOC1
Species	Human
GeneID (Human)	846
GeneID	846
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 211-310 of human CaSR (NP_000379.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P41180
Protein Weight	121kDa
Shipping	Ice bag