


Pe Rabbit Anti-Human Cd100/Sema4D Monoclonal Antibody For Flow Cytometry - Q92854

Item Number: CM0010400



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD100/SEMA4D, validated for flow cytometry. Recognizes human Semaphorin-4D, a cell surface receptor involved in immune and neuronal signaling. Ideal for immunology and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD100/SEMA4D mAb
Remarks/Alias	A8; GR3; BB18; CD100; COLL4; SEMA4; coll-4; C9orf164; M-sema-G
Species	Human
GeneID (Human)	10507
GeneID	10507
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 22-734 of human Semaphorin-4D/CD100, NP_006369.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q92854
Protein Weight	82kDa/96kDa
Shipping	Ice bag