

Pe Rabbit Anti-Human Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P10721

Item Number: CM0010377



Introduction

PE-conjugated rabbit monoclonal antibody against human CD117/c-Kit, validated for flow cytometry. Enables detection of the SCF receptor in hematopoiesis, mast cell, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD117/c-Kit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC
Species	Human
GeneID	3815
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-520 of human CD117/c-Kit (NP_000213.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag