


Pe Rabbit Anti-Human Cd11A Mab For Fc - P20701

Item Number: CM0009568



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD11a (ITGAL, LFA-1) for flow cytometry (FC). Validated for human samples, with molecular weights 118-134 kDa. Supports immune cell adhesion and leukocyte trafficking research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD11a mAb
Remarks/Alias	CD11A; LFA-1; LFA1A
Species	Human
GeneID (Human)	3683
GeneID	3683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a (NP_002200.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20701
Protein Weight	118kDa / 129kDa / 134kDa
Shipping	Ice bag