


Pe Rabbit Anti-Human Cd123 Mab - P26951

Item Number: CM0004072



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD123 (IL-3 receptor alpha). Validated for flow cytometry, this antibody detects native CD123 on hematopoietic cells, supporting immunophenotyping and leukemia diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD123 mAb
Remarks/Alias	IL3R; CD123; IL3RX; IL3RY; IL3RAY; hIL-3Ra; IL-3R-alpha
Species	Human
GeneID (Human)	3563
GeneID	3563
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-378 of human CD123 (NP_002174.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26951
Protein Weight	43 kDa
Shipping	Ice bag