

Pe Rabbit Anti-Human Cd124/IL-4Rα Monoclonal Antibody For Flow Cytometry - P24394

Item Number: CM0006254



Introduction

This PE-conjugated rabbit monoclonal antibody recognizes Human CD124 (IL-4Rα) and is validated for flow cytometry. It is suitable for studies on Th2 differentiation, IgE production, and allergic inflammation responses.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human CD124/IL-4Rα mAb
Remarks/Alias	CD124; IL4RA; IL-4RA
Species	Human
GeneID (Human)	3566
GeneID	3566
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-232 of human CD124/IL-4Rα (NP_000409.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24394
Protein Weight	26kDa/88kDa/90kDa
Shipping	Ice bag