

Pe Rabbit Anti-Human Cd131/Csf2Rb Monoclonal Antibody For Flow Cytometry - P32927

Item Number: CM0007803



Introduction

PE-conjugated rabbit monoclonal antibody against human CD131 (CSF2RB), validated for flow cytometry. Targets the common beta subunit of IL-3/IL-5/GM-CSF receptors. Suitable for hematopoietic and immune cell research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD131 mAb
Remarks/Alias	CD131; IL3RB; IL5RB; SMDP5; CDw131; betaGMR
Species	Human
Gene ID (Human)	1439
Gene ID	1439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-443 of human CD131 (NP_000386.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32927
Protein Weight	97kDa
Shipping	Ice bag