


Pe Rabbit Anti-Human Cd140B/Pdgfrβ Monoclonal Antibody For Flow Cytometry - P09619

Item Number: CM0010135



Introduction

This PE-conjugated rabbit monoclonal antibody targets human CD140b/PDGFRβ (PDGFRB) and is validated for direct flow cytometry. It is suitable for studying receptor signaling in angiogenesis, cancer biology, and vascular development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD140b/PDGFRβ mAb
Remarks/Alias	PDGFRB; CD140B; IBGC4; IMF1; JTK12; KOGS; PDGFR; PDGFR-1; PDGFR1; PENTT; platelet-derived growth factor receptor beta
Species	Human
Gene ID (Human)	5159
Immunogen	Recombinant fusion protein corresponding to amino acids 33-530 of human CD140b/PDGFRB (NP_002600.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09619
Protein Weight	37kDa/123kDa
Shipping	Ice bag