

Pe Rabbit Anti-Human Cd141/Thrombomodulin Mab - P07204

Item Number: CM0010804



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD141/Thrombomodulin, validated for flow cytometry. Recombinant immunogen (aa 19-515). Suitable for endothelial cell, coagulation, and dendritic cell (BDCA-3) profiling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD141/Thrombomodulin mAb
Remarks/Alias	TM; THRM; AHUS6; BDCA3; CD141; BDCA-3; THPH12
Species	Human
GeneID (Human)	7056
GeneID	7056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-515 of human CD141/Thrombomodulin (NP_000352.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07204
Protein Weight	60kDa
Shipping	Ice bag