

Pe Rabbit Anti-Human Cd146 Monoclonal Antibody For Fc - P43121

Item Number: CM0009402



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD146 (MCAM/MUC18). Designed for flow cytometry, it reacts with human samples. Ideal for studying cell adhesion, endothelial biology, and melanoma metastasis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD146 mAb
Remarks/Alias	MCAM; CD146; MUC18; cell surface glycoprotein MUC18
Species	Human
GeneID (Human)	4162
GeneID	4162
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-559 of human CD146 (NP_006491.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43121
Protein Weight	57kDa/71kDa
Shipping	Ice bag