

Pe Rabbit Anti-Human Cd152/Ctla-4 Monoclonal Antibody For Intracellular Flow Cytometry - P16410

Item Number: CM0005271



Introduction

PE-conjugated rabbit monoclonal antibody against human CD152/CTLA-4, an inhibitory T-cell receptor. Designed for intracellular flow cytometry (FC) applications. Reacts specifically with human CD152/CTLA-4.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD152/CTLA-4 mAb
Remarks/Alias	CD; GSE; GRD4; ALPS5; CD152; CTLA-4; IDDM12; CELIAC3
Species	Human
GeneID (Human)	1493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-162 of human CD152/CTLA-4 (NP_005205.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P16410
Protein Weight	6kDa/8kDa/19kDa/24kDa
Shipping	Ice bag