

## Pe Rabbit Anti-Human Cd153 Mab - P32971

Item Number: CM0005932

### Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD153 (CD30 ligand), validated for flow cytometry. Suitable for T-cell proliferation and lymphoma research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human CD153 mAb                                                                                    |
| Remarks/Alias    | CD153; CD30L; CD30LG; TNLG3A                                                                                      |
| Species          | Human                                                                                                             |
| GeneID (human)   | 944                                                                                                               |
| GeneID           | 944                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 63-234 of human CD153 (NP_001235.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | P32971                                                                                                            |
| Protein Weight   | 26kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |