


Pe Rabbit Anti-Human Cd156C/Adam10 Monoclonal Antibody For Flow Cytometry - O14672

Item Number: CM0011676



Introduction

PE-conjugated rabbit monoclonal antibody targeting human ADAM10 (CD156c), validated for flow cytometry. ADAM10 is a transmembrane metalloprotease involved in ectodomain shedding of proteins like TNF, APP, and Notch, with roles in inflammation, neurodegeneration, and cancer.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD156c/ADAM10 mAb
Remarks/Alias	RAK; kuz; AD10; AD18; MADM; CD156c; CDw156; HsT18717
Species	Human
GeneID (Human)	102
GeneID	102
Immunogen	Recombinant protein corresponding to a sequence within amino acids 214-672 of human ADAM10 (NP_001101.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14672
Protein Weight	84kDa
Shipping	Ice bag