

Pe Rabbit Anti-Human Cd158A/Kir2DL1 Monoclonal Antibody For Flow Cytometry - P43626

Item Number: CM0009690



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD158a/KIR2DL1 (NK cell inhibitory receptor). Validated for flow cytometry (FC), cross-reacts with human. Suitable for NK cell immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD158a/KIR2DL1 mAb
Remarks/Alias	NKAT; NKAT1; p58.1; CD158A; KIR221; NKAT-1; KIR-K64; KIR2DL3
Species	Human
Gene ID (Human)	3802
Gene ID	3802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human CD158a/KIR2DL1 (NP_055033.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43626
Protein Weight	39kDa
Shipping	Ice bag