

Pe Rabbit Anti-Human Cd161 Mab For Flow Cytometry - Q12918

Item Number: CM0010731



Introduction

PE-conjugated rabbit monoclonal antibody against human CD161 (KLRB1/NKR-P1A) for flow cytometry. Reacts with human. Useful for studying NK cell inhibitory receptor function and immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD161 mAb
Remarks/Alias	NKR; CD161; CLEC5B; NKR-P1; NKR-P1A; NKR-P1A; hNKR-P1A
Species	Human
Gene ID	3820
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-225 of human CD161 (NP_002249.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12918
Protein Weight	25kDa
Shipping	Ice bag