


Pe Rabbit Anti-Human Cd163 Mab For Fc - Q86Vb7

Item Number: CM0009156



Introduction

PE-conjugated rabbit monoclonal antibody against human CD163 (M130, hemoglobin scavenger receptor) for flow cytometry. Validated for human samples, this antibody is ideal for macrophage and monocyte studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD163 mAb
Remarks/Alias	M130; MM130; SCAR11
Species	Human
GeneID (human)	9332
GeneID	9332
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 42-1045 of human IL-10RA/CD210 (NP_004235.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q86VB7
Protein Weight	125 kDa
Shipping	Ice bag