


Pe Rabbit Anti-Human Cd166/Alcam Monoclonal Antibody For Flow Cytometry - Q13740

Item Number: CM0008174



Introduction

PE-conjugated rabbit monoclonal antibody specifically recognizing human CD166/ALCAM (SWISS: Q13740). Developed for flow cytometry, it supports studies on cell adhesion, immune synapse formation, stem cell homing, and tumor metastasis.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human CD166/ALCAM mAb
Remarks/Alias	MEMD; CD166
Species	Human
Gene ID (Human)	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-526 of human CD166/ALCAM (NP_001618.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13740
Protein Weight	15kDa/33kDa/63kDa/65kDa
Shipping	Ice bag