

Pe Rabbit Anti-Human Cd167A/Ddr1 Monoclonal Antibody For Fc - Q08345

Item Number: CM0010222



Introduction

PE-conjugated rabbit monoclonal antibody against human CD167a/DDR1, validated for flow cytometry (FC). Epithelial discoidin domain receptor 1, a collagen-activated tyrosine kinase, implicated in cancer cell invasion and tissue remodeling.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD167a/DDR1 mAb
Remarks/Alias	CAK; DDR; NEP; HGK2; PTK3; RTK6; TRKE; CD167; EDDR1; MCK10; NTRK4; PTK3A
Species	Human
GeneID (Human)	780
GeneID	780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-416 of human CD167a/DDR1 (NP_054699.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q08345
Protein Weight	101kDa
Shipping	Ice bag