


## Pe Rabbit Anti-Human Cd170/Siglec-5 Mab - O15389

Item Number: CM0011982



### Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD170/Siglec-5, a sialic-acid-binding lectin expressed on myeloid cells. Validated for flow cytometry, suitable for immunology and inflammation research.

[Learn More](#)

| Parameter        | Value                                   |
|------------------|-----------------------------------------|
| Product Name     | PE Rabbit anti-Human CD170/Siglec-5 mAb |
| Remarks/Alias    | CD170, OBBP2, CD33L2, OB-BP2, SIGLEC-5  |
| Species          | Human                                   |
| GeneID (Human)   | 8778                                    |
| GeneID           | 8778                                    |
| Immunogen        | Recombinant protein                     |
| Source           | Rabbit                                  |
| Category         | Monoclonal Antibodies                   |
| Application      | FC                                      |
| Cross Reactivity | Human                                   |
| SWISS            | O15389                                  |
| Protein Weight   | 61kDa                                   |
| Shipping         | Ice bag                                 |