


Pe Rabbit Anti-Human Cd172A/Sirpα Monoclonal Antibody For Flow Cytometry - P78324

Item Number: CM0011184



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD172a/SIRPα (P78324), validated for flow cytometry. Suitable for studying phagocytosis, immune regulation, and synaptic signaling mechanisms.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD172a/SIRPα mAb
Remarks/Alias	BIT; MFR; P84; SIRP; MYD-1; SHPS1; CD172A; PTPNS1
Species	Human
GeneID (Human)	140885
GeneID	140885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-370 of human CD172a/SIRPα (NP_542970.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P78324
Protein Weight	55 kDa
Shipping	Ice bag