

Pe Rabbit Anti-Human Cd177 Monoclonal Antibody For Fc - Q8N6Q3

Item Number: CM0005371

Introduction

PE-conjugated rabbit monoclonal antibody against human CD177 (NB1/PRV-1) for flow cytometry. Designed to detect CD177 on neutrophil subsets; cross-reacts with human. Suitable for immunophenotyping and research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD177 mAb
Remarks/Alias	NB1; PRV1; HNA2A; PRV-1; HNA-2a; NB1 GP
Species	Human
GeneID (human)	57126
GeneID	57126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-407 of human CD177 (NP_065139.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8N6Q3
Protein Weight	15kDa/27kDa/46kDa
Shipping	Ice bag