

Pe Rabbit Anti-Human Cd185/Cxcr5 Monoclonal Antibody For Flow Cytometry - P32302

Item Number: CM0008972

Introduction

PE-conjugated rabbit monoclonal antibody against human CD185/CXCR5 (BLR1), validated for flow cytometry. Targets the B-cell chemoattractant receptor involved in B-cell migration and Burkitt lymphoma. Swiss-Prot: P32302.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD185/CXCR5 Monoclonal Antibody
Remarks/Alias	BLR1; CD185; MDR15
Species	Human
GeneID (Human)	643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-47aa of human CD185/CXCR5 (NP_033981.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32302
Protein Weight	36kDa/42kDa
Shipping	Ice bag