


Pe Rabbit Anti-Human Cd197/Ccr7 Monoclonal Antibody For Flow Cytometry - P32248

Item Number: CM0012539



Introduction

PE-conjugated rabbit monoclonal antibody against human CD197/CCR7 (BLR2, EBI1) for flow cytometry. Suitable for studying lymphocyte trafficking, homing, and Epstein-Barr virus-induced receptor expression. UniProt P32248.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD197/CCR7 mAb
Remarks/Alias	BLR2; EBI1; CCR-7; CD197; CDw197; CMKBR7; CC-CKR-7
Species	Human
GeneID (Human)	1236
GeneID	1236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-59 of human CD197/CCR7 (NP_001829.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32248
Protein Weight	43kDa
Shipping	Ice bag