


Pe Rabbit Anti-Human Cd1B Monoclonal Antibody For Flow Cytometry - P29016

Item Number: CM0012083



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD1b (SWISS: P29016), validated for flow cytometry, ideal for studies on lipid antigen presentation and T-cell immunity.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD1B mAb
Remarks/Alias	R1; CD1; CD1A
Species	Human
GeneID (Human)	910
GeneID	910
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human CD1B (NP_001755.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P29016
Protein Weight	37kDa/30kDa
Shipping	Ice bag