


Pe Rabbit Anti-Human Cd1C Monoclonal Antibody For Flow Cytometry - P29017

Item Number: CM0008126



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD1c (R7, BDCA1), a lipid antigen-presenting molecule. Validated for flow cytometry, this 38 kDa antibody is ideal for studies on adaptive immunity and glycolipid presentation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD1c mAb
Remarks/Alias	R7; CD1; CD1A; BDCA1
Species	Human
GeneID (Human)	911
GeneID	911
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-298 of human CD1c (NP_001756.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P29017
Protein Weight	38kDa
Shipping	Ice bag