

Pe Rabbit Anti-Human Cd2 Monoclonal Antibody For Flow Cytometry - P06729

Item Number: CM0011736



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD2 (T11, LFA-2, SRBC) specifically for flow cytometry. Ideal for T-cell and NK cell identification in immunological research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD2 mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
GeneID (Human)	914
GeneID	914
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag