


Pe Rabbit Anti-Human Cd200R Monoclonal Antibody For Flow Cytometry - Q8Td46

Item Number: CM0008739



Introduction

PE-conjugated rabbit monoclonal antibody against human CD200R (OX2 receptor) validated for flow cytometry. Ideal for immunophenotyping and exploring immune checkpoint pathways in inflammation and oncology research. Molecular weight: 37 kDa; species reactivity: human.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD200R mAb
Remarks/Alias	OX2R; MOX2R; CD200R; HCRTR2
Species	Human
Gene ID (Human)	131450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-240 of human CD200R (NP_740750.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TD46
Protein Weight	37kDa
Shipping	Ice bag