


Pe Rabbit Anti-Human Cd205/Dec-205 Monoclonal Antibody For Fc - O60449

Item Number: CM0010928



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD205/DEC-205 (LY75) for flow cytometry. Detects the endocytic receptor on dendritic cells and B cells. High specificity, validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD205/DEC-205 mAb
Remarks/Alias	CD205; LY-75; CLEC13B; DEC-205; GP200-MR6
Species	Human
GeneID (Human)	4065
GeneID	4065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-501 of human CD205/DEC-205 (NP_002340.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O60449
Protein Weight	198kDa
Shipping	Ice bag