

Pe Rabbit Anti-Human Cd207/Langerin Mab For Fc - Q9Uj71

Item Number: CM0005359



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD207/Langerin for flow cytometry. Specific for human, ideal for studying Langerhans cell biology and antigen presentation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD207/Langerin mAb
Remarks/Alias	CLEC4K
Species	Human
GeneID	50489
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-328 of human Langerin/CD207 (NP_056532.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UJ71
Protein Weight	37kDa
Shipping	Ice bag