


Pe Rabbit Anti-Human Cd213A1/IL-13R α 1 Monoclonal Antibody For Flow Cytometry - P78552

Item Number: CM0009987



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD213a1/IL-13R α 1 (CT19) for flow cytometry. Reacts with human. 49 kDa. Ideal for studying interleukin-13 signaling in immune cells.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD213a1/IL-13R α 1 mAb
Remarks/Alias	NR4; CT19; CD213A1; IL-13Ra
Species	Human
GeneID (Human)	3597
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 29-343 of human CD213a1/IL-13R α 1 (NP_001551.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P78552
Protein Weight	49kDa
Shipping	Ice bag