

Pe Rabbit Anti-Human Cd220 Monoclonal Antibody For Flow Cytometry - P06213

Item Number: CM0008234



Introduction

PE-conjugated rabbit monoclonal antibody against human CD220 (Insulin receptor) for flow cytometry. Recognizes human CD220 with a molecular weight of 156 kDa. Validated for FC research use.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD220 mAb
Remarks/Alias	HHF5; CD220
Species	Human
GeneID (Human)	3643
GeneID	3643
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-944 of human CD220 (NP_000199.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06213
Protein Weight	156 kDa
Shipping	Ice bag