

Pe Rabbit Anti-Human Cd221A/Igf-1Rα Mab - P08069

Item Number: CM0011157



Introduction

PE-conjugated rabbit anti-human CD221/IGF-1R monoclonal antibody validated for flow cytometry. Suitable for detecting IGF-1R expression in human cells, supporting cancer and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD221α/IGF-1Rα mAb
Remarks/Alias	IGFR; CD221; IGFR; JTK13
Species	Human
GeneID (Human)	3480
GeneID	3480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-932 of human CD221α/IGF-1Rα (NP_000866.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08069
Protein Weight	155kDa
Shipping	Ice bag