


## Pe Rabbit Anti-Human Cd226/Dnam-1 Monoclonal Antibody For Fc - Q15762

Item Number: CM0009341



### Introduction

PE-conjugated rabbit anti-human CD226 (DNAM-1) monoclonal antibody for flow cytometry. This antibody targets CD226, a costimulatory receptor on T and NK cells that boosts cytotoxicity against tumors and infected cells upon binding PVR/CD155 and Nectin2/CD112. Ideal for immune cell profiling.

[Learn More](#)

| Parameter        | Value                                                                                                                                          |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human CD226/DNAM-1 mAb                                                                                                          |
| Remarks/Alias    | PTA1; DNAM1; DNAM-1; TLISA1                                                                                                                    |
| Species          | Human                                                                                                                                          |
| GeneID (Human)   | 10666                                                                                                                                          |
| GeneID           | 10666                                                                                                                                          |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 19-247 of human CD226/DNAM-1 (NP_006557.2) |
| Source           | Rabbit                                                                                                                                         |
| Category         | Monoclonal Antibodies                                                                                                                          |
| Application      | FC                                                                                                                                             |
| Cross Reactivity | Human                                                                                                                                          |
| SWISS            | Q15762                                                                                                                                         |
| Protein Weight   | 39kDa                                                                                                                                          |
| Shipping         | Ice bag                                                                                                                                        |