

Pe Rabbit Anti-Human Cd23 Monoclonal Antibody For Fc - P06734

Item Number: CM0010129



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD23 (FcεRII) for flow cytometry. Detects low-affinity IgE receptor involved in B cell regulation, with cross-reactivity to human and molecular weight of 36 kDa.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD23 mAb
Remarks/Alias	CD23; FCE2; CD23A; IGEBF; CLEC4J; FcεRII; BLAST-2
Species	Human
Gene ID (Human)	2208
Gene ID	2208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-321 of human CD23 (NP_001993.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06734
Protein Weight	36kDa
Shipping	Ice bag