

Pe Rabbit Anti-Human Cd24 Monoclonal Antibody For Flow Cytometry - P25063

Item Number: CM0004190

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD24, suitable for flow cytometry. CD24 is a GPI-anchored glycoprotein involved in cell differentiation and immune modulation. Validated cross-reactivity with human; ideal for B-cell studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD24 mAb
Remarks/Alias	CD24A; CD24
Species	Human
GeneID (Human)	100133941
GeneID	100133941
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-80 of human CD24 (NP_037362.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25063
Protein Weight	8kDa; 13kDa
Shipping	Ice bag