

Pe Rabbit Anti-Human Cd244/2B4 Monoclonal Antibody For Fc - Q9Bzw8

Item Number: CM0008626

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD244/2B4 (SLAMF4) for flow cytometry applications. Recognizes human CD244, a key NK cell receptor, suitable for immune cell profiling and activation studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD244/2B4 mAb
Remarks/Alias	CD244; 2B4; NAIL; NKR2B4; Nmrk; SLAMF4; CD244 molecule
Species	Human
GenelD (human)	51744
GenelD	51744
Immunogen	Recombinant protein corresponding to a sequence within amino acids 22-221 of human CD244/2B4 (NP_001160135.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BZW8
Protein Weight	30kDa/36kDa/41kDa
Shipping	Ice bag