


Pe Rabbit Anti-Human Cd262/Dr5/Trailr2 Monoclonal Antibody For Flow Cytometry - O14763

Item Number: CM0010188



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD262/DR5/TRAILR2, optimized for flow cytometry. Detects the TRAIL receptor in human samples. Ideal for apoptosis and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD262/DR5/TRAILR2 mAb
Remarks/Alias	DR5; CD262; KILLER; TRICK2; TRICKB; ZTNFR9; TRAILR2; TRICK2A; TRICK2B; TRAIL-R2; KILLER/DR5
Species	Human
GeneID (human)	8795
GeneID	8795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 56-182 of human CD262/DR5/TRAILR2 (NP_003833.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14763
Protein Weight	48kDa
Shipping	Ice bag