

## Pe Rabbit Anti-Human Cd263/Trail-R3 Monoclonal Antibody For Fc - O14798

Item Number: CM0004557

### Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD263/TRAIL-R3 (decoy receptor 1). Designed for flow cytometry (FC) to detect TRAIL-R3 on cell surfaces. Facilitates apoptosis and cancer research. - O14798

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human CD263/TRAIL-R3 mAb                                                                                    |
| Remarks/Alias    | LIT; DCR1; TRID; CD263; TRAILR3; TRAIL-R3; DCR1-TNFR                                                                       |
| Species          | Human                                                                                                                      |
| GenelD (human)   | 8794                                                                                                                       |
| GenelD           | 8794                                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-235 of human CD263/TRAIL-R3(NP_003832.3). |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | FC                                                                                                                         |
| Cross Reactivity | Human                                                                                                                      |
| SWISS            | O14798                                                                                                                     |
| Protein Weight   | 27kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |