

## Pe Rabbit Anti-Human Cd270/Hvem Monoclonal Antibody For Fc - Q92956

Item Number: CM0011683



### Introduction

PE-conjugated rabbit monoclonal antibody against human CD270/HVEM (TNFRSF14), validated for flow cytometry. Suitable for analyzing T cell costimulation, immune checkpoint signaling, and herpesvirus entry in human samples.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human CD270/HVEM mAb                                                                                    |
| Remarks/Alias    | TNFRSF14; ATAR; CD270; HVEA; HVEM; LIGHTR; TR2; TNF receptor superfamily member 14                                     |
| Species          | Human                                                                                                                  |
| Gene ID (Human)  | 8764                                                                                                                   |
| Gene ID          | 8764                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 39-162 of human CD270/HVEM (NP_003811.2) |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | FC                                                                                                                     |
| Cross Reactivity | Human                                                                                                                  |
| SWISS            | Q92956                                                                                                                 |
| Protein Weight   | 21kDa/30kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                                |