

Pe Rabbit Anti-Human Cd275/B7-H2/Icosl Monoclonal Antibody For Fc - O75144

Item Number: CM0005824



Introduction

PE-conjugated rabbit monoclonal antibody against human CD275/B7-H2/ICOSL, specifically designed for flow cytometry (FC). Ideal for studying T-cell costimulation, B-cell activation, and neutrophil transmigration in immune research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD275/B7-H2/ICOSL mAb
Remarks/Alias	B7h; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
Species	Human
Gene ID (Human)	23308
Gene ID	23308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human CD275/B7-H2/ICOSL (NP_056074.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75144
Protein Weight	20kDa/33kDa/34kDa
Shipping	Ice bag